

MCM2 Rabbit mAb


ABZ4946

Product Overview


Name	MCM2 Rabbit mAb
Accession(Primary)	P49736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:200; ICC/IF-1:200-1:1000; FC-1:100-1:200

Target information (P49736)

Gene Name	MCM2
About Gene Symbol	



Western blot analysis of MCM2 in HeLa lysates using MCM2 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse colon using MCM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.