

ABCG1 Rabbit mAb


ABZ4947

Product Overview

Name	ABCG1 Rabbit mAb
Accession(Primary)	P45844
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 76 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:100

Target information (P45844)

Gene Name	ABCG1
About Gene Symbol	



Western blot analysis of ABCG1 in HUVEC lysates
using ABCG1 antibody.