

CD134 Receptor Rabbit mAb

ABZ4948

Product Overview

Name	CD134 Receptor Rabbit mAb
Accession(Primary)	P43489
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 29 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:20

Target information (P43489)

Gene Name TNFRSF4**About Gene Symbol**

