

PI3 Kinase p110 beta Rabbit mAb


ABZ4950

Product Overview

Name	PI3 Kinase p110 beta Rabbit mAb
Accession(Primary)	P42338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 123 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (P42338)

Gene Name	PIK3CB
About Gene Symbol	PIK3CB (Phosphatidylinositol-4,5-Bisphosphate 3-Kinase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely s



Western blot analysis of PI3 Kinase p110 beta in Jurkat lysates using PI3 Kinase p110 beta antibody.