

CD79B Rabbit mAb


ABZ4953

Product Overview

Name	CD79B Rabbit mAb
Accession(Primary)	P40259
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 30-50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:20-1:50; FC-1:20-1:50

Target information (P40259)

Gene Name	CD79B
About Gene Symbol	



Western blot analysis of CD79B in Ramos, rat spleen lysates using CD79B antibody.