


Thrombopoietin Receptor Rabbit mAb

ABZ4954

Product Overview

Name	Thrombopoietin Receptor Rabbit mAb
Accession(Primary)	P40238
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P40238)

Gene Name MPL**About Gene Symbol**

Western blot analysis of TPOR in (1)A431 lysates; (2)NIH/3T3 lysates using Thrombopoietin Receptor antibody.

