

MMP12 Rabbit mAb


ABZ4955

Product Overview

Name	MMP12 Rabbit mAb
Accession(Primary)	P39900
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54,45,22 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P39900)

Gene Name	MMP12
About Gene Symbol	



Western blot analysis of MMP12 in WI38 lysates
using MMP12 antibody.