

p21 Rabbit mAb


ABZ4956

Product Overview


Name	p21 Rabbit mAb
Accession(Primary)	P38936
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 18 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:500-1:1000; FC-1:50-1:100; IP-1:10-1:100

Target information (P38936)

Gene Name	CDKN1A
About Gene Symbol	



Western blot analysis of p21 Waf1/Cip1 in MCF-7 lysates using p21 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using p21 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.