

ZEB1 Rabbit mAb


ABZ4957

Product Overview

Name	ZEB1 Rabbit mAb
Accession(Primary)	P37275
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 124 kDa; Observed MW: 200 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (P37275)

Gene Name	ZEB1
About Gene Symbol	ZEB1 (Zinc Finger E-Box Binding Homeobox 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ZEB1/TCF8 in Jurkat, HT-1080, Ramos lysates using ZEB1 antibody.