

ATP6V1E1 Rabbit mAb


ABZ4958

Product Overview


Name	ATP6V1E1 Rabbit mAb
Accession(Primary)	P36543
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1; 50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P36543)


Gene Name	ATP6V1E1
About Gene Symbol	



Western blot analysis of ATP6V1E1 in 3T3, HeLa, 293 lysates using ATP6V1E1 antibody.



Immunocytochemistry analysis of ATP6V1E1(green) in HeLa using ATP6V1E1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ATP6V1E1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.