

VEGF Receptor 2 Rabbit mAb


ABZ4959

Product Overview

Name	VEGF Receptor 2 Rabbit mAb
Accession(Primary)	P35968
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 152 kDa; Observed MW: 230,210 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P35968)

Gene Name	KDR
About Gene Symbol	



Western blot analysis of VEGF Receptor 2 in Human placenta lysates using VEGF Receptor 2 antibody.