

STIP1/STI1 Rabbit mAb


ABZ4962

Product Overview

Name	STIP1/STI1 Rabbit mAb
Accession(Primary)	P31948
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; FC-1:20-1:50

Target information (P31948)

Gene Name	STIP1
About Gene Symbol	



Western blot analysis of STIP1 in Jurkat lysates
using STIP1 antibody.