

RRM2 Rabbit mAb


ABZ4965

Product Overview

Name	RRM2 Rabbit mAb
Accession(Primary)	P31350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200; IP-1:10-1:100

Target information (P31350)

Gene Name	RRM2
About Gene Symbol	



Western blot analysis of RRM2 in HeLa lysates
using RRM2 antibody.