

PRKAR2B Rabbit mAb


ABZ4966

Product Overview

Name	PRKAR2B Rabbit mAb
Accession(Primary)	P31323
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:20-1:50; IP-1:20-1:50

Target information (P31323)

Gene Name	PRKAR2B
About Gene Symbol	



Western blot analysis of PKA 2 beta in Human fetal brain lysates using PKA RII beta antibody.