


# Angiotensin II Type 1 Receptor Rabbit mAb

ABZ4967

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Angiotensin II Type 1 Receptor Rabbit mAb  |
| <b>Accession(Primary)</b>   | P30556                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 41 kDa; Observed MW: 41 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                           |

## Target information (P30556)

**Gene Name** AGTR1**About Gene Symbol**

Western blot analysis of AGTR1 in (1) HeLa lysates; (2) PC-12 lysates using Angiotensin II Type 1 Receptor antibody.

