

Aquaporin 1 Rabbit mAb


ABZ4969

Product Overview

Name	Aquaporin 1 Rabbit mAb
Accession(Primary)	P29972
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:200-1:2000

Target information (P29972)

Gene Name	AQP1
About Gene Symbol	AQP1 (Aquaporin 1 (Colton Blood Group)) is a formin actin nucleator involved in adaptive immunity. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of AQP1 in Human fetal kidney lysates using Aquaporin 1 antibody.