


SHP1 Rabbit mAb

ABZ4970


Product Overview

Name	SHP1 Rabbit mAb
Accession(Primary)	P29350
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P29350)

Gene Name PTPN6**About Gene Symbol**

Western blot analysis of SHP1 in Jurkat lysates using SHP1 antibody.



Immunohistochemistry analysis of paraffin-embedded rat spleen using SHP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.