

PTBP1 Rabbit mAb


ABZ4972

Product Overview

Name	PTBP1 Rabbit mAb
Accession(Primary)	P26599
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:200-1:2000

Target information (P26599)

Gene Name	PTBP1
About Gene Symbol	PTBP1 (Polypyrimidine Tract Binding Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PTBP1 in Daudi lysates
using PTBP1 antibody.