

HMGB2 Rabbit mAb


ABZ4973

Product Overview

Name	HMGB2 Rabbit mAb
Accession(Primary)	P26583
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500

Target information (P26583)

Gene Name	HMGB2
About Gene Symbol	HMGB2 (High Mobility Group Box 2) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HMGB2 in (1) HeLa
lysates; (2) PC-12 lysates using HMGB2 antibody.