

YY1 Rabbit mAb


ABZ4975

Product Overview


Name	YY1 Rabbit mAb
Accession(Primary)	P25490
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:100; IF-1:20-1:50; FC-1:20-1:50

Target information (P25490)


Gene Name	YY1
About Gene Symbol	YY1 (YY1 Transcription Factor) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of YY1 in (1) HeLa lysates;
(2) Daudi lysates using YY1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human bladder using YY1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of YY1 in HeLa using YY1 antibody.