

Stathmin 1 Rabbit mAb


ABZ4980

Product Overview


Name	Stathmin 1 Rabbit mAb
Accession(Primary)	P16949
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (P16949)

Gene Name	STMN1
About Gene Symbol	



Western blot analysis of Stathmin 1 in C6, HeLa, SH-SY5Y lysates using Stathmin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Stathmin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.