

BNP Rabbit mAb


ABZ4981

Product Overview

Name	BNP Rabbit mAb
Accession(Primary)	P16860
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 15 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:50-1:100

Target information (P16860)

Gene Name	NPPB
About Gene Symbol	



Western blot analysis of BNP in Human plasma
lysates using BNP antibody.