

S100 beta Rabbit mAb


ABZ4982

Product Overview


Name	S100 beta Rabbit mAb
Accession(Primary)	P04271
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100

Target information (P04271)


Gene Name	S100B
About Gene Symbol	



Western blot analysis of S100B in A375 lysates using S100B antibody.



Immunohistochemistry analysis of paraffin-embedded rat brain using S100B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of S100B in A375 using S100B antibody.

