

ER alpha Rabbit mAb


ABZ4986

Product Overview


Name	ER alpha Rabbit mAb
Accession(Primary)	P03372
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:100-1:200; ICC/IF-1:100-1:200; FC-1:20

Target information (P03372)

Gene Name	ESR1
About Gene Symbol	



Western blot analysis of Estrogen Receptor α in MCF-7 lysates using Estrogen Receptor alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Estrogen Receptor alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.