

Human IgM Rabbit mAb


ABZ4989

Product Overview


Name	Human IgM Rabbit mAb
Accession(Primary)	P01871
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 49 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; IP-1:20-1:50

Target information (P01871)

Gene Name	IGHM
About Gene Symbol	



Western blot analysis of Human IgM in Human plasma lysates using Human IgM antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Human IgM antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.