


Human IgG3 Rabbit mAb

ABZ4990

Product Overview

Name	Human IgG3 Rabbit mAb
Accession(Primary)	P01860
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (P01860)

Gene Name IGHG3**About Gene Symbol**

Western blot analysis of Human IgG3 in Human plasma lysates using Human IgG3 antibody.

