

CD8 alpha Rabbit mAb


ABZ4991

Product Overview


Name	CD8 alpha Rabbit mAb
Accession(Primary)	P01732
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500; IHC-1:50-1:100

Target information (P01732)


Gene Name	CD8A
About Gene Symbol	




Western blot analysis of CD8 in Jurkat lysates using CD8 alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using CD8 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of CD8 alpha in Jurkat using CD8 antibody.



Western blot analysis of CD8 alpha in Raji lysates using CD8 antibody.