

TGF alpha Rabbit mAb


ABZ4992

Product Overview

Name	TGF alpha Rabbit mAb
Accession(Primary)	P01135
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P01135)

Gene Name	TGFA
About Gene Symbol	



Western blot analysis of TGF alpha in (1) HepG2 lysates; (2) Human fetal liver lysates using TGF alpha antibody.