

MTCO2 Rabbit mAb


ABZ4994

Product Overview


Name	MTCO2 Rabbit mAb
Accession(Primary)	P00403
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 26 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (P00403)

Gene Name	MT-CO2
About Gene Symbol	



Western blot analysis of MTCO2 in K562 lysates using MTCO2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using MTCO2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.