

Flotillin 1 Rabbit mAb


ABZ5001

Product Overview


Name	Flotillin 1 Rabbit mAb
Accession(Primary)	O75955
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:10-1:100; IP-1:10-1:100

Target information (O75955)

Gene Name	FLOT1
About Gene Symbol	FLOT1 (Flotillin 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Flotillin 1 in (1) HeLa lysates; (2) K562 lysates using Flotillin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human uterus using Flotillin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.