

MKLP1 Rabbit mAb


ABZ5003

Product Overview

Name	MKLP1 Rabbit mAb
Accession(Primary)	Q02241
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 110 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:200-1:500

Target information (Q02241)

Gene Name	KIF23
About Gene Symbol	KIF23 (Kinesin Family Member 23) is a kinesin motor protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MKLP1 in HeLa lysates using MKLP1 antibody.