

Caveolin-2 Rabbit mAb


ABZ5005

Product Overview


Name	Caveolin-2 Rabbit mAb
Accession(Primary)	P51636
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:100-1:200; IF-1:100-1:200; FC-1:200-1:2000

Target information (P51636)

Gene Name	CAV2
About Gene Symbol	



Western blot analysis of Caveolin 2 in HeLa lysates using Caveolin 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using Caveolin2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.