

ERK2 Rabbit mAb


ABZ5006

Product Overview


Name	ERK2 Rabbit mAb
Accession(Primary)	P28482
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:200-1:2000; FC-1:200-1:1000

Target information (P28482)

Gene Name	MAPK1
About Gene Symbol	MAPK1 (Mitogen-Activated Protein Kinase 1) is a receptor tyrosine kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of ERK2 in (1) A431 lysates;
(2) HeLa lysates using ERK2 antibody.



Immunohistochemistry analysis of paraffin-
embedded Human gastric adenocarcinoma using
ERK2 antibody. High-pressure and temperature
Sodium Citrate pH 6.0 was used for antigen
retrieval.