

PDGFR beta Rabbit mAb


ABZ5007

Product Overview


Name	PDGFR beta Rabbit mAb
Accession(Primary)	P09619
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 124 kDa; Observed MW: 170 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20

Target information (P09619)


Gene Name	PDGFRB
About Gene Symbol	



Western blot analysis of PDGF Receptor beta in SH-SY5Y lysates using PDGFR beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human uterus using PDGF Receptor beta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of PDGFR beta in 3T3 using PDGF Receptor beta antibody.

