

KDEL Rabbit mAb


ABZ5008

Product Overview


Name	KDEL Rabbit mAb
Accession(Primary)	P24390
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 78KD kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (P24390)

Gene Name	KDELRL1
About Gene Symbol	



Western blot analysis of KDEL in HeLa lysates using KDEL antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using KDEL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.