


CSNK2A1 alpha Rabbit mAb

ABZ5010


Product Overview

Name	CSNK2A1 alpha Rabbit mAb
Accession(Primary)	P68400
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:20-1:50; IHC-1:50-1:100

Target information (P68400)

Gene Name CSNK2A1**About Gene Symbol**

Western blot analysis of CKII alpha in HeLa lysates using CKII alpha antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using CKII alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.