

c-Met Rabbit mAb


ABZ5013

Product Overview


Name	c-Met Rabbit mAb
Accession(Primary)	P08581
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 156 kDa; Observed MW: 156 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; IF-1:50-1:200

Target information (P08581)


Gene Name	MET
About Gene Symbol	



Western blot analysis of c-Met in 293 lysates using c-Met antibody.



Immunohistochemistry analysis of paraffin-embedded Human stomach using Met (c-Met) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of c-Met in HT-29 using Met (c-Met) antibody.

