

Thrombomodulin Rabbit mAb


ABZ5016

Product Overview


Name	Thrombomodulin Rabbit mAb
Accession(Primary)	P07204
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:1000; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:10-1:100

Target information (P07204)

Gene Name	THBD
About Gene Symbol	



Western blot analysis of Thrombomodulin in Human placenta lysates using Thrombomodulin antibody.



Immunohistochemistry analysis of paraffin-embedded Human bladder using Thrombomodulin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.