

ACADM/MCAD Rabbit mAb


ABZ5017

Product Overview

Name	ACADM/MCAD Rabbit mAb
Accession(Primary)	P11310
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:200; FC-1:50-1:100; IP-1:10-1:100

Target information (P11310)

Gene Name	ACADM
About Gene Symbol	ACADM (Acyl-CoA Dehydrogenase Medium Chain 2 3 5) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mech



Western blot analysis of ACADM in (1) HeLa lysates; (2) K562 lysates using ACADM antibody.