

Nibrin Rabbit mAb

ABZ5018


Product Overview

Name	Nibrin Rabbit mAb
Accession(Primary)	O60934
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200


Target information (O60934)

Gene Name NBN

About Gene Symbol



Western blot analysis of p95/NBS1 in HeLa lysates using Nibrin antibody.



Immunohistochemistry analysis of paraffin-embedded rat heart using p95/NBS1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.