

Bad Rabbit mAb


ABZ5019

Product Overview


Name	Bad Rabbit mAb
Accession(Primary)	Q92934
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20; FC-1:50-1:100

Target information (Q92934)


Gene Name	BAD
About Gene Symbol	



Western blot analysis of BAD in Hela lysates using Bad antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Bad antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Bad in Hela using Bad antibody.

