

DARPP32 Rabbit mAb


ABZ5020

Product Overview

Name	DARPP32 Rabbit mAb
Accession(Primary)	Q9UD71
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; IP-1:20

Target information (Q9UD71)

Gene Name	PPP1R1B
About Gene Symbol	



Western blot analysis of DARPP32 in rat brain
lysates using DARPP32 antibody.