

PKC beta 2 Rabbit mAb


ABZ5021

Product Overview

Name	PKC beta 2 Rabbit mAb
Accession(Primary)	P05771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; FC-1:20-1:50

Target information (P05771)

Gene Name	PRKCB
About Gene Symbol	



Western blot analysis of PKC beta 2 in K562
lysates using PKC beta 2 antibody.