


Human IgG4 Rabbit mAb

ABZ5023


Product Overview

Name	Human IgG4 Rabbit mAb
Accession(Primary)	P01861
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:500-1:5000

Target information (P01861)

Gene Name IGHG4**About Gene Symbol**

Western blot analysis of IgG4 in Human spleen lysates using Human IgG4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human stomach carcinoma using IgG4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.