

Interferon alpha/beta Receptor 1 Rabbit mAb


ABZ5024

Product Overview

Name	Interferon alpha/beta Receptor 1 Rabbit mAb
Accession(Primary)	P17181
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 64 kDa; Observed MW: 90,130 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:50-1:100; IP-1:20-1:50

Target information (P17181)

Gene Name	IFNAR1
About Gene Symbol	



Western blot analysis of IFNAR1 in K562 lysates
using Interferon alpha/beta Receptor 1 antibody.