

Integrin beta 1 Rabbit mAb


ABZ5025

Product Overview


Name	Integrin beta 1 Rabbit mAb
Accession(Primary)	P05556
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000

Target information (P05556)

Gene Name	ITGB1
About Gene Symbol	ITGB1 (Integrin Subunit Beta 1) is an integrin cell adhesion receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Integrin beta 1 in A431 lysates using Integrin beta 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Integrin beta 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.