


Histone H2B (yeast) Rabbit mAb

ABZ5026

Product Overview

Name	Histone H2B (yeast) Rabbit mAb
Accession(Primary)	P02293
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Saccharomyces cerevisiae
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P02293)

Gene Name HTB1**About Gene Symbol**

Western blot analysis of Histone H2B in
Saccharomyces cerevisiae lysates using Histone
H2B (yeast) antibody.

