

IRF1 Rabbit mAb


ABZ5028

Product Overview


Name	IRF1 Rabbit mAb
Accession(Primary)	P10914
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:100-1:200; IP-1:50-1:100

Target information (P10914)


Gene Name	IRF1
About Gene Symbol	



Western blot analysis of IRF1 in Jurkat lysates using IRF1 antibody.



Immunohistochemistry analysis of paraffin-embedded rat kidney using IRF1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of IRF1 in HeLa using IRF1 antibody.