

Histone H2A Rabbit mAb


ABZ5030

Product Overview


Name	Histone H2A Rabbit mAb
Accession(Primary)	P04908
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:2000; ICC/IF-1:100-1:200; FC-1:200-1:1000

Target information (P04908)

Gene Name	H2AC4
About Gene Symbol	



Western blot analysis of Histone H2A.Z
in(1)Neuro2a lysates;(2)HeLa lysates using Histone
H2A antibody.



Immunohistochemistry analysis of paraffin-
embedded Human kidney using Histone H2A.Z
antibody.High-pressure and temperature Sodium
Citrate pH 6.0 was used for antigen retrieval.