

HIF1 alpha Rabbit mAb


ABZ5031

Product Overview

Name	HIF1 alpha Rabbit mAb
Accession(Primary)	Q16665
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 93 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; IHC:1:50-1:100; ICC/IF-1:200-1:500; FC-1:500-1:2000; IP-1:50-1:200

Target information (Q16665)

Gene Name	HIF1A
About Gene Symbol	



Western blot analysis of HIF1alpha in HL-60
lysates using HIF1 alpha antibody.